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DOCUMENT 3

Trypsin (from bovine pancreas) EC 3.4.21.4Catalysed reaction :

Hydrolyses proteins, peptides, amides and esters specifically at the carboxyl groups of the basic amino acids L -arginine and L- lysine.

Characteristics :

Molecular weight : ca. 23 000

Isoelectric point : 10.5

pH-optimum : 8.0

Inhibitors : Natural inhibitors from pancreas, soybean, limabean and egg with ; organophosphorus compounds such as diisopropyl fluorophosphate ; several amines and amidines ; calcium.

Stability :

Store dry at 4°C ; no considerable loss of activity over a period of 12 months.

DOCUMENT 4

E G T A

Ethylene glycol-bis (β -aminoethylether) N, N, N', N' -tetra acetic acid.
A chelating agent useful for the determination of calcium in the presence of magnesium.

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


Glycerol

UV-method

for the determination of glycerol in foodstuffs and other materials

Cat.No. 148270


Test-Combination
for 3 x 10 determinations

Principle


In the reaction catalyzed by glycerokinase (GK), glycerol is phosphorylated by adenosine-5'-triphosphate (ATP) to glycerol-3-phosphate (1).

The adenosine-5'-diphosphate (ADP) formed in the above reaction is reconverted by phosphoenolpyruvate (PEP) with the aid of pyruvate kinase (PK) into ATP with the formation of pyruvate (2).

In the presence of the enzyme lactate dehydrogenase (LDH) pyruvate is reduced to lactate by reduced nicotinamide-adenine dinucleotide (NADH) with the oxidation of NADH to NAD (3).

The amount of NADH consumed in the above reaction is stoichiometric with the amount of glycerol. NADH is determined by means of its absorption at 334, 340 or 365 nm.

Each Test-Combination Contains

- Three bottles 1 with 2 g coenzyme/buffer mixture each, consisting of
glycylglycine buffer — pH 7.4
NADH — 7 mg
ATP — 22 mg
PEP — 11 mg
magnesium sulphate
stabilizers
- Bottle 2 with 0.4 ml enzyme suspension, consisting of
pyruvate kinase — ca. 240 U
lactate dehydrogenase — ca. 220 U
- Bottle 3 with 0.4 ml enzyme suspension
glycerokinase — ca. 34 U

Preparation of Solutions for 10 Determinations

- Dissolve contents of one bottle 1 with 11 ml redist. water. Before use let the solution stand for about 10 min at room temperature.
- Use contents of bottle 2 undiluted.
- Use contents of bottle 3 undiluted.

Stability of Solutions

Solution 1 is stable for 4 days at +4°C.

Suspensions of bottles 2 and 3 are stable for 1 year at +4°C.

Procedure

Wavelength¹: 340 nm, Hg 365 nm or Hg 334 nm

Glass cuvette²: 1 cm light path

Temperature: 20–25°C

Final volume: 3.02 ml


Read against air (without cuvette in the light path) or against water.

Sample solution: 3–50 µg glycerol/cuvette³ (in 0.1–2.0 ml sample volume).

¹ The absorption maximum of NADH is at 340 nm. With spectrophotometers, measurements are made at the absorption maximum; when employing spectrum line photometers equipped with a mercury vapour lamp, measurements are made at a wavelength of 365 nm or 334 nm.

² Disposable cuvettes may be used instead of glass cuvettes, if desired.

³ See instructions for performance of assay and sample preparation.

Recommendations to methods and standardized procedures see references.

pipette into cuvettes	blank	sample
solution 1	1.00 ml	1.00 ml
redist. water	2.00 ml	1.90 ml
sample solution	—	0.10 ml
suspension 2	0.01 ml	0.01 ml
mix, after completion of the reaction (ca. 5–7 min) read absorbances of the solutions (A ₁). Start reaction by addition of		
suspension 3	0.01 ml	0.01 ml
Mix, wait for the completion of the reaction (ca. 5–10 min) and read the absorbances of the solutions (A ₂). If the reaction has not stopped after 15 min, continue to read the absorbances at 2 min intervals until the absorbance decreases constantly.		

If the absorbance at A₂ decreases constantly, extrapolate the absorbances to the time of the addition of suspension 3.

Calculate the absorbance differences (A₁–A₂) for the blank and the sample. Subtract the absorbance difference of the blank from the absorbance difference of the sample.

$$\Delta A = \Delta A_S - \Delta A_B$$

Calculation

According to the general formula for calculating the concentration the equation is:

$$c = \frac{V \times MW}{\epsilon \times d \times v \times 1000} \times \Delta A \text{ [g/l]}$$

where V = final volume [ml]

v = sample volume [ml]

MW = molecular weight of the substance to be assayed

d = light path [cm]

ε = absorption coefficient of NADH at:

$$340 \text{ nm} = 6.3 \text{ [l} \times \text{mmol}^{-1} \times \text{cm}^{-1}]$$

$$\text{Hg } 365 \text{ nm} = 3.4 \text{ [l} \times \text{mmol}^{-1} \times \text{cm}^{-1}]$$

$$\text{Hg } 334 \text{ nm} = 6.18 \text{ [l} \times \text{mmol}^{-1} \times \text{cm}^{-1}]$$

It follows for glycerol

$$c = \frac{3.02 \times 92}{\epsilon \times 1 \times 0.1 \times 1000} \times \Delta A = 2.781 \times \frac{\Delta A}{\epsilon} \text{ [g glycerol/l sample solution]}$$

If the sample was diluted during preparation, the result must be multiplied by the dilution factor F.